




**RENTMED**

**SPANDAN**

**ECG AT HOME**



## Revolutionizing Cardiac Care with Portable ECG Technology

- **Where Cable Meets Intelligence**  
Spandan Pro simplifies ECGs by combining traditional leads with smart technology.
- **Direct Smartphone Connectivity**  
Spandan Pro's direct smartphone connection makes the ECG procedure convenient.
- **Eliminating Bulky ECG Machines**  
Focusing on the portability and user-friendliness of Spandan Pro, making traditional ECG machines obsolete.
- **Portability for On-the-Go Diagnosis**  
Enabling cardiologists and physicians to conduct high-quality ECGs anywhere, from clinics to remote patient visits.

Highly detailed reports with accurate data

### 12 Lead ECG Report

Date: 16 Feb 2024, 10:11 AM  
 Name: Vignesh Kumar  
 Age: 34 years  
 Gender: Male  
 Height: 175 cm  
 Weight: 68 kg  
 Heart Rate: 60 bpm

#### ECG Parameters

Parameter	Observed Value	Standard Range
PR Interval	155 ms	100 ms - 200 ms
QRS Interval	98 ms	60 ms - 120 ms
QT Interval	379 ms	300 ms - 450 ms
QTc Interval	457 ms	300 ms - 450 ms
Heart Rate	60 bpm	60 bpm - 120 bpm

#### Result Details

Overall evaluation of test:  
 Normal ECG  
 Interpretation details:  
 Normal  
 Anomalies in test reports:  
 Negative  
 Suggestive action:  
 Follow healthy life style

#### Risk Meter



### ECG Parameters

Parameter	Observed Value	Standard Range
PR Interval	155 ms	100 ms - 200 ms
QRS Interval	98 ms	60 ms - 120 ms
QT Interval	379 ms	300 ms - 450 ms
QTc Interval	457 ms	300 ms - 450 ms
Heart Rate	60 bpm	60 bpm - 120 bpm


#### Risk Meter



### 12 Lead ECG Report

Date: 16 Feb 2024, 10:11 AM  
 Name: Vignesh Kumar  
 Age: 34 years  
 Gender: Male  
 Height: 175 cm  
 Weight: 68 kg  
 Heart Rate: 60 bpm

#### 12 Lead ECG Report



## 10 Auto-Switching Chest Leads



No need to move the leads



Leads record Sequential ECG



3D Manual for easy positioning




Fast Interpretations



Leads with high tensile strength




Consistent signal transfer






# CONTACT US



[www.rentmed.co.in](http://www.rentmed.co.in)




+91 9520990890

# OFFICE ADDRESS

**Address: 19/20, Prakash Nagar, Chowk  
Number 1, Govindgarh Road, Dehradun,  
Uttarakhand 248001, India**

## 2 Types Of ECG Interpretation

### Interpretation by Spandan's AI



Instant  
Reports



Highly  
Detailed



Easy to  
Understand

MAKE IN INDIA

### Interpretation by Spandan's Experts



**Step-1** Take ECG & generate report

**Step-2** WhatsApp us the reports

**Step-3** Get manual interpretation!

Disclaimer: This service is available on a subscription.



# Product Description

## Cutting Edge Technology

Spandan features bio-potential sensors and Intelli-ECG tech for 12-lead, 7-lead, and lead II ECG tests, 5-minute HRV test, Live ECG monitoring, and Framingham's Heart Risk Calculator. FDA-cleared with 99.7% trace accuracy.

## Comprehensive Detection


Detects AFib, Bradycardia, Tachycardia, Atrial Flutter, AV-Block, Ventricular Tachycardia, and other heart abnormalities like QRS variations, Myocardial Infarctions, P wave absence, and T wave changes.

## Record and Share On The Go

Spandan ECG device records preliminary tests within 10 seconds anytime, anywhere. One can easily record, save, share unlimited ECG reports via your smartphone. Compatible with all android devices.

# Goldberger ECG System

Spandan Pro includes 10 leads and features the Goldberger ECG System, which is an alternative to the standard ECG system.





## Cutting-Edge Technology in Cardiac Diagnostics

